

Reusable Biopsy Forceps

/Pinces à biopsie à usage multiple /Pinzas de biopsia reutilizables



PERFORMANCE[®]

Highly polished stainless steel jaws laser welded

/Mors en acier hautement poli et soudure laser /Mandíbulas en acero inoxidable pulidas y soldadas con laser

Fully polished spiral without welding avoiding damage to the endoscope

/Spirale entièrement polie sans soudure afin de limiter les risques d'endommager l'endoscope/Espiral enteramente pulida sin soldaduras evitando dañar el endoscopio

Large cups permitting a better biopsy sample

/Grande cuillère permettant l'obtention d'une meilleure biopsie/Copa de grande tamaño permitiendo una mejor muestra de biopsia

Internal closing mechanism without hook avoiding damage to the endoscope

/Mécanisme interne de fermeture sans crochet limitant les risques d'endommager l'endoscope/Sistema de cierre interno sin ganchos evitando dañar el endoscopio

Serial number for complete traceability

/Traçabilité complète par numéro de lot unique /Número de serie permitiendo su seguimiento



Pediatric /Pédiatrique /Pediátrica



Bronco /Bronchique /Bronco



Gastro /Gastrique /Gastro



Colono /Colique /Colono

Exclusive ergonomic color identification handle

/Poignée ergonomique avec code couleur exclusif
/Mango ergonómico con color para identificación

Reusable Biopsy Forceps

/Pincas à biopsie à usage multiple /Pinzas de biopsia reutilizables

PERFORMANCE



Reference Réf. / Ref.	Length (cm) Long. / Longitud	Diam. (mm) Diam. / Diám.	Cup opening (mm) Ouverture des mors / Apertura de la copa	Cup length (mm) Longueur des mors / Longitud de la copa	Description Description / Descripción
GF120	120	1.8	5.9	2.8	Bronco
GF160	160	2.3	7.4	3.4	Gastro
GF160P	160	1.8	5.9	2.8	Gastro/Pediatric
GF230	230	2.3	7.4	3.4	Colono
GF230P	230	1.8	5.9	3.3	Colono/Pediatric
GF23034	230	3.4	8.3	4.6	Colono



Reference Réf. / Ref.	Length (cm) Long. / Longitud	Diam. (mm) Diam. / Diám.	Cup opening (mm) Ouverture des mors / Apertura de la copa	Cup length (mm) Longueur des mors / Longitud de la copa	Description Description / Descripción
GF160S	160	2.3	7.1	3.4	Gastro With Spike
GF230S	230	2.3	7.1	3.4	Colono With Spike
GF160PS	160	1.8	5.9	2.8	Gastro, Pediatric with Spike



Cup opening
/Ouverture des mors
/Apertura de la copa



Cup length
/Longueur des mors
/Longitud de la copa